Instruction Manual for LDB Series

LDB Series Despatch Ovens are bench ovens to 204°C (400°F) with forced convection airflow.

MODEL	<u>VOLTS</u>	HEATER WATTS	<u>AMPS</u>	<u>HZ</u>	<u>PHASE</u>
LDB 1-17 LDB 1-24 LDB 1-43 LDB 1-69 LDB 2-18 LDB 2-27	120** 120** 120** 120** 240	1,200 1,200 1,600 2,400 3,600 4,800	11.6 11.6 15.0 21.6 16.6 21.6	50/60 50/60 50/60 50/60 50/60	1 1 1 1 1

^{**} A 240V conversion kit is available.

Notice

Users of this equipment must comply with operating procedures and training of operation personnel as required by the Occupational Safety and Health Act (OSHA) of 1970, Section 6 and relevant safety standards, as well as other safety rules and regulations of state and local governments. Refer to the relevant safety standards in OSHA and National Fire Protection Association (NFPA), section 86 of 1990.

Caution

Setup and maintenance of the equipment should be performed by qualified personnel who are experienced in handling all facets of this type of system. Improper setup and operation of this equipment could cause an explosion that may result in equipment damage, personal injury or possible death.

Dear Customer,

Thank you for choosing Despatch Industries. We appreciate the opportunity to work with you and to meet your heat processing needs. We believe that you have selected the finest equipment available in the heat processing industry.

At Despatch, our service does not end after the purchase and delivery of our equipment. For this reason we have created the Service Products Division within Despatch. The Service Products Division features our Response Center for customer service. The Response Center will direct and track your service call to ensure satisfaction.

Whenever you need service or replacement parts, contact the Response Center at 1-800-473-7373: FAX 612-781-5353.

Thank you for choosing Despatch.

Sincerely,

Despatch Industries

PREFACE

This manual is your guide to the Despatch oven. It is organized to give you the information you need quickly and easily.

The INTRODUCTION section provides an overview of the Despatch oven.

The THEORY OF OPERATION section details the function and operation of assemblies and subassemblies on the Despatch oven.

The INSTRUCTIONS section provides directions on unpacking, installing, operating and maintaining the Despatch oven.

The APPENDIX section contains Special Instructions for operating the control instrument, a Troubleshooting Table, a list of Accessories and a Warranty.

An efficient way to learn about the oven would be to read the manual while working with the corresponding oven control system. This will give you practical hands-on experience with information in the manual and the oven.

Before operating the equipment, be sure you understand all of the technical information contained in this manual. Information skipped, not understood or misunderstood could create the possibility of operating the equipment in an unsafe manner. This can cause damage to the oven or personnel or reduce the efficiency of the equipment.

NOTE:
Read the entire
INTRODUCTION and
THEORY OF OPERATION
before installing the oven.

WARNING:

Failure to heed warnings in this instruction manual and on the oven could result in personal injury, property damage or death.

Table of Contents

PREFACE	
Table of Contents	
List of Figures	iv
List of Tables	iv
INTRODUCTION	
Special Features	. 1
Specifications	. 2
Dimensions	
Power	
Temperature	. 4
Capacities	. 4
THEORY OF OPERATION	
CONTROL Instrument	
HI-LIMIT Instrument	_
Product HI-LIMIT Instrument	_
Oven HI-LIMIT Instrument	
Oven Theory	. 9
INSTRUCTIONS	. 11
Unpacking and Inspection	. 11
Setup	. 12
HI-LIMIT Instrument Adjustment	. 13
Operating	. 14
Loading the Oven	. 14
Pre-Startup Checklist	. 13
Startup	, 10 10
Maintenance	, 10 10
Checklist	10
Tests	19
Parts	19
CONTROL Instrument	
HI-LIMIT Instrument	21
Heater Unit	22
Fan Motor	23
A DDFN IDIV	25
APPENDIX	. 25 25
Special Instructions	. 25 26
Calibration	. 20
CONTROL	. 30
Calibration Instructions	30

Table of Contents (Cont.)

HI-LIMIT	. 31
Calibration Instructions	31
Accessories	. 32
Troubleshooting	. 33
Warranty	. 35
Despatch Classic [™] Service Warranty Program	. 35
Immediate Service Response	35
Field Service Network	35
Best Service Protection in the Industry!	36
Classic [™] Series Warranty	. 37
Parts and Material	37
Labor and Expenses	37
Exclusions/Limitations of Liability	37
Emergency Service	38
Non-Compliance	- 38
Despatch Customer Service	. 39
Procedures and Customer Responsibilities	. 39
Attachment A - Sustained Service Support	. 40
Drawings	. 41
LDB1-17 120V	. 41
LDB1-17 240V	. 42
List of Material	43
LDB1-24 120V	. 44
LDB1-24 240V	. 45
List of Material	46
LDB1-43 120V	. 47
LDB1-43 240V	. 48
List of Material	49
LDB1-69 120V	. 50
LDB1-69 240V	. 51
List of Material	52
LDB2-18 240V	. 53
List of Material	54
LDB2-27 240V	. 55
List of Material	56

List of Figures

Figure 1 illustrates the forced circulating airflow in the Despatch
LDB ovens
Figure 2 illustrates the CONTROL instrument
Figure 3 illustrates the HI-LIMIT instrument.
Figure 4 illustrates the control panel on the LDB Series oven.
Figure 5 illustrates the connections to the CONTROL
instrument
Figure 6 illustrates the connections to the HI-LIMIT instrument.
Figure 7 illustrates the control instrument
Figure 8 illustrates the Classic ™4 Plus 1 Service Program 30
Liet of Tobles
List of Tables
Table 1 Dimensions
Table 2 Power Requirements
Table 3 Temperature Specifications
Table 4 Capacities
Table 5 CONTROL Instrument Features
Table 6 Motor Leads to Terminal Strip2
iable of Moror Leads to Terrinal Suit
Table 7 Parameter Program Mode Outline

INTRODUCTION

The INTRODUCTION section provides an overview of the Despatch LDB Series Ovens. The LDB Series Ovens have the most effective heat distribution and the fastest processing time of any lab oven its size. Air is discharged from the floor plate of the oven and is distributed through the chamber.

Special Features

The sturdy construction and high grade insulation of the Despatch LDB Series ovens contribute to excellent high temperature performance. Other special features include:

- Dual functioning (Proportional or ON-OFF) digital CONTROL instrument to control temperature fluctuations.
- Digital CONTROL and manual reset HI-LIMIT instrument to protect the chamber workload as well as the oven itself.
- Unique Despatch design to combine forced circulated air with a perforated stainless steel floor plate for the ultimate in temperature performance.
- Welded double wall construction and fiberglass insulation to reduce heat loss. Silicone rubber gaskets further minimize heat leakage.
- Rapid response heater with a five year warranty.
- Scratch-resistant Silver-Clay® baked enamel exterior and stainless steel interior for easy cleaning.
- Space-saving, stackable design.

Specifications

Dimensions

Table 1 Dimensions

Model	690000000000000000000000000000000000000	mber n (cm		Capacity feet		verall Si in (cm)	ze	Shelves Provided on	Maximum Number of	Chamber
1VIOGE1	W	D	Н	(liters)	W	D	Н	Shelf Centers	Shelf Positions	Doors
LDB 1-18	18 (46)	14 (36)	12 (30)	1.75 (4 9.5)	21 (53)	19 (48)	24 (61)	2 on 2"	4	1
LDB 1-24	24 (61)	14 (36)	12 (30)	2.3 (65)	27 (69)	19 (48)	24 (61)	2 on 2"	4	1
LDB 1-43	24 (61)	14 (36)	22 (56)	4.3 (122)	27 (69)	19 (48)	35 (89)	2 on 2"	9	1
LDB 1-69	30 (76)	18 (46)	22 (56)	6.9 (195)	34 (86)	23 (58)	36 (89)	2 on 2"	9	2
LDB 2-18	37 (94)	24 (61)	35 (89)	18 (510)	40 (102)	29 (74)	49 (122)	2 on 2"	15	2
LDB 2-27	37 (94)	37 (94)	35 (89)	27.7 (785)	40 (102)	42 (107)	49 (122)	2 on 2"	15	2

Power

Line voltages may vary in some geographical locations. If your line voltage is much lower than the oven voltage rating, warm up time will be longer and motors may overload or run hot. If your line voltage is higher than name plate rating, the motor may run hot and draw excessive amps.

If the line voltage varies more than 10% from the oven voltage rating, some electrical components such as relays, temperature controls, etc. may operate erratically.

Table 2 Power Requirements

Model	Volts	Amps	Hertz	Phase	Heater (KW)	Cord and Plug
LDB 1-17	120 240*	11.6 5.8	50/60	1	1.2	Included, 15 Amp
LDB 1-24	120 240*	11.6 5.8	50/60	1	1.2	Included, 15 Amp
LDB 1-43	120 240*	15.0 7.5	50/60	1	1.6	Included, 20 Amp
LDB 1-69	120 240*	21.6 10.8	50/60	1	2.4	None, Hardwired
LDB 2-18	240	16.6	50/60	1	3.6	None, Hardwired
LDB 2-27	240	21.6	50/60	1	4.8	None, Hardwired

The LDB 1-69, LDB 2-18 and LDB 2-27 must be hardwired to the electric supply using 10AWG or larger wires suitable for at least 75°C (167°F).

Ovens designed for 240 volts (see name plate on oven) will operate satisfactorily on a minimum of 208 volts, but with a 25% reduction in heater power. If your power characteristics are lower, contact Despatch Industries.

^{*} Order 240 conversion kit for 240 on LDB 1-17, LDB 1-24, LDB 1-43 and LDB 1-69.

Temperature

Table 3 Temperature Specifications

				LI)B		
Model		1-17	1-24	1-43	1-69	2-18	2-27
Time to Temperature (approximate minutes with no load)	25°C - 100°C 25°C - 150°C 25°C - 204°C	6 11 22	8 14 29	10 19 31	9 17 26	8 16 26	9 17 28
Recovery Time Door Open 1 Min. (approximation minutes with no load)	100°C nate 150°C 204°C	1 3 5	1 3 5	1 3 5	2 3 5	2 3 5	2 3 5
Operating Range w/20°C Ambient			40	50°C -	204°C		

Capacities

Table 4 Capacities

	LDB						
Model	1-17	1-24	1-43	1-69	2-18	2-27	
Maximum load capacity (lbs)	75	75	75	100	150	200	
Maximum shelf load (lbs)	25	30	25	20	50	50	
Exhaust Capacity, CFM	2½	2½	2½	6	14	14	
Fan CFM H.P.	150 1/25	150 1/25	150 1/25	150 1/25	300 2/25 (2)1/25	300 2/25 (2)1/25	

THEORY OF OPERATION

The THEORY OF OPERATION section details the function and operation of assemblies and subassemblies on the Despatch LDB Series Ovens.

The Despatch LDB Series Ovens have the most effective heat distribution system and the fastest processing time of any lab oven its class. These ovens are especially useful for testing, preheating, sterilizing, drying, aging and curing as well as other production applications.

The Despatch LDB Series Ovens incorporate forced circulating airflow with precision digital control to deliver fast processing. The overall result is efficient productivity under strenuous conditions. The LDB Series Ovens are precise and practical.




Figure 1 illustrates the forced circulating airflow in the Despatch LDB ovens.

The unique Despatch design, moves convected heat through a perforated stainless steel floor plate. The air is circulated with a high volume fan. The Despatch LDB Series Ovens employ higher volume fans than any competitive model. The chamber can be densely loaded without interfering with the process. Air temperature near the sensing element in the chamber is within 1°C of the process temperature appearing on the digital display. For your convenience the fresh air intake is fixed. The exhaust rate is regulated by dampers on the top of the unit.

CONTROL Instrument

The oven is equipped with a microprocessor based digital control instrument.




Figure 2 illustrates the CONTROL instrument.

Table 5 CONTROL Instrument Features

Features	Description
Main Display	Displays the actual oven temperature or displays the setpoint when the set key is pressed.
Program Key and ewtr910 Key	Change operating mode parameters from factory preset values. Allow access to program mode when pressed simultaneously with set key.
Set Key	Views the setpoint or, in conjunction with the ▲ and ▼ keys, changes the setpoint.
Down Key (▼)	Decreases a setpoint or mode parameter.
Up Key (*)	Increases a setpoint or mode parameter.
LED Output Indicator	Lights when the control is calling for heat.
LED SV Indicator	Lights when the setpoint value is displayed.

CONTROL Instrument (Cont.)

The Despatch CONTROL instrument is a dual functioning proportional or ON-OFF controller. Due to the physical characteristics of the oven, the CONTROL has been configured as an ON-OFF controller and set to its optimum operating values. Initially the CONTROL will allow the heater to operate at full power. However, as the actual oven temperature reaches the setpoint, the ON-OFF CONTROL will cycle the heater on and off, minimizing process temperature fluctuations.

HI-LIMIT Instrument

The oven is equipped with a HI-LIMIT instrument. The purpose of the HI-LIMIT instrument is to provide a protective measure for the product and/or the oven itself. If the setting on the HI-LIMIT is exceeded, the heating process will discontinue, thus protecting the product and/or the oven.

Set the HI-LIMIT instrument to a temperature 10°C - 14°C higher than the CONTROL instrument setpoint or a temperature that should not be exceeded in the process. If the setting on the HI-LIMIT instrument is exceeded the heater will shut down. The HI-LIMIT instrument must be manually reset by pushing the reset button on the HI-LIMIT instrument.



Figure 3 illustrates the HI-LIMIT instrument.

Product HI-LIMIT Instrument

If the product being processed has a critical high temperature limit, the HI-LIMIT instrument should be used as a product HI-LIMIT instrument. The HI-LIMIT instrument should be set to a temperature somewhat below the temperature at which the product could be damaged. Use the CONTROL instrument or a pyrometer to determine the product high limit setting.

Oven HI-LIMIT Instrument

If the product does not have a critical high temperature limit, the HI-LIMIT can be used as an oven HI-LIMIT instrument. An oven HI-LIMIT instrument protects oven equipment. Since the HI-LIMIT instrument does not show the temperature, it can be accurately set only after the oven is in operation. Until then, the HI-LIMIT should be set at the maximum position so all preliminary testing and adjusting can be done.

Oven Theory

The oven has an efficient forced circulating oven to 204°C (400°F). A forced circulating oven relies on a circulating motor to move air through the chamber. A forced circulating chamber is much more efficient and uniform than a gravity convected oven. However, the amount of airflow in a circulating oven is low compared to ovens with a forced recirculating airflow system. Also, it takes a finite amount of time for the oven to soak in at the desired setpoint. The time that it takes the unit to soak in at setpoint is related to such parameters as chamber area, load mass and the ability to absorb heat and exhaust rate. The unique Despatch design and CONTROL action compensates for these facts.

The oven uses an indicating microprocessor based digital control that displays the actual chamber temperature at the sensing point. Despatch strategically locates the CONTROL's temperature sensor to optimize the control action for the entire chamber for various load conditions. The CONTROL display will fluctuate a few degrees around the setpoint, but the overall chamber temperature will remain very stable. The underlying reason for this is that the display is showing temperature fluctuations at the temperature sensor location, not the overall chamber temperature. The strategic location of the sensor compensates for delays in heat convection and enhances the performance and temperature control of the oven. The oven has been engineered to have an overall result of quality productivity where fast processing and versatility are critical.

INSTRUCTIONS

The INSTRUCTIONS section provides directions on unpacking, installation, operation and maintenance of the Despatch LDB Series Ovens.

Unpacking and Inspection

Remove all packing materials and thoroughly inspect the oven for damage of any kind that could have occurred during shipment.

- See whether the carton and plastic cover sheet inside carton are still in good condition.
- Look at all outside surfaces and corners of the oven for scratches and dents.
- Check the oven controls and indicators for normal movement, bent shafts, cracks, chips or missing parts such as knobs and lenses.
- · Check the door and latch for smooth operation.

If there is damage, and it could have happened during shipment follow these instructions.

- 1. Contact the shipper immediately and file a written damage claim.
- Contact Despatch Industries to report your findings and to order replacement parts for those that were damaged or missing.
- 3. Please send a copy of your filed damage claims to Despatch.

Unpacking and Inspection (Cont.)

Next, check to make sure you have received all the required materials. Your shipment should include:

- One (1) Despatch oven,
- One (1) Instruction manual,
- One (1) Warranty card,
- Two (2) Shelves
- One (1) Package containing four rubber feet

If any of these items are missing from the packaged contents, contact Despatch Industries to have the appropriate materials forwarded to you.

Finally, to protect the warranty on your new LDB Series Oven, complete the warranty card and mail it to Despatch within 15 days after receipt of the equipment.

Setup

- 1. Remove adhesive backing sheet from the rubber feet.
- 2. Attach rubber feet to the bottom corners of the oven.
- 3. Place oven on a bench top or an optional cabinet base.

The oven must have a minimum of two (2) inches clearance in the rear to provide proper ventilation. The oven may be placed next to another cabinet, or next to another oven, with three (3) inch clearance (the doors will still open).

Make sure oven is level and plumb; this will assure proper heat distribution and operation of all mechanical components.

- 4. Identify correct power source indicated on the specification plate.
- 5. Plug or hardwire oven directly to the electric supply.

WARNING:

All grounding and safety equipment must be in compliance with applicable codes, ordinances and accepted safe practices.

HI-LIMIT Instrument Adjustment

Before putting the oven into production, adjust the oven HI-LIMIT instrument as follows.

- 1. Press the **POWER** switch to ON.
- 2. Set the CONTROL instrument at 14°C (25°F) above the desired operating temperature or the desired high limit temperature.
- 3. Push the black button to reset the HI-LIMIT instrument and operate the oven until the CONTROL instrument is regulating.
- 4. Allow the oven to soak at this temperature for 30 minutes.
- 5. Carefully adjust the HI-LIMIT instrument downward until it trips. The **HEATER ON** LED will come on.
- 6. Reset the CONTROL instrument at the desired operating temperature.

It will be necessary to reset the HI-LIMIT instrument whenever it has tripped. The HI-LIMIT instrument may be reset by first allowing the oven chamber to cool slightly (or by tuning the HI-LIMIT instrument thermostat up several degrees) and pushing the black reset button.

Operating

Users and operators of this oven must comply with operating procedures and training of operating personnel as required by the Occupational Safety and Health Act (OSHA) of 1970, Section 5 and relevant safety standards, and other safety rules and regulations of state and local governments. Refer to the relevant safety standards in OSHA and National Fire Protection Association (NFPA), Section 86 of 1990.

WARNING:

Do not use oven in wet, corrosive or explosive atmospheres unless this oven is specifically designed for a special atmosphere.

Loading the Oven

Despatch Industries cannot be responsible for either the process or process temperature used, or for the quality of the product being processed. It is the responsibility of the purchaser and operator to see that the product undergoing processing in a Despatch oven is adequately protected from damage.

Carefully following the instructions in this manual will help the purchaser and operator in fulfilling that responsibility.

When loading the oven avoid spills of anything onto the heater elements or onto the floor of the oven. Do not place the load on the oven floor plate. This may cause the load to heat unevenly and the weight may cause shorting out of the heater elements. Use the shelves provided.

The two shelves are designed to be pulled out about half way without tipping. The support capacity of the shelves is listed in the Capacities Table in the Specifications section in this manual. Do not overload the shelves.

Distribute the workload evenly so that airflow is not restricted. Do not overfill your oven. The workload should not take up more than two-thirds of any dimension of the inside cavity.

WARNING:

Never operate oven at a temperature in excess of the maximum operating temperature of 204°C (400°F).

Pre-Startup Checklist

- ✓ <u>Know the system.</u> Read this manual carefully. Make use of its instructions and explanations. The know how of safe, continuous, satisfactory, trouble-free operation depends primarily on the degree of your understanding of the system and of your willingness to keep all parts in proper operating condition.
- ✓ <u>Check line voltage</u>. Voltage must correspond to nameplate requirements of motors and controls. Refer to the section on power connections in the INTRODUCTION of this manual.
- ✓ Fresh air and exhaust. Do not be careless about restrictions in and around the fresh air and exhaust openings. Under no condition permit them to become so filled with dirt that they appreciably reduce the air quantity. Refer to the Set-up instructions in this manual.
- ✓ <u>Ventilation</u> There is an exhaust opening in the top of the unit.

The exhaust vent(s) may have to be closed to reach the maximum temperature of 204°C, especially if operating on 208 volts. They may need to be opened to operate properly at the lower range of the oven's design

✓ Helpful hints

- For drying ovens, open vent to prevent buildup of moisture.
- For sample heating, close the vent when no ventilation is required.

NOTE:

The vents may have to be adjusted to achieve maximum performance at various operating temperatures.

Startup

For fastest oven heat-up time, close the fresh-air vent. After the desired temperature is reached, the vent may be adjusted as needed.

- 1. Start oven.
 - a. Open oven door.
 - b. Press **POWER ON** pushbutton. You will hear the circulating fan start.
 - c. Shut oven door.
 - d. Set vent to desired opening. The vents may have to be adjusted to achieve maximum performance at various operating temperatures.
 - e. Check that the amber LED above the **POWER** switch is on.




Figure 4 illustrates the control panel on the LDB Series oven.

- 2. Enter setpoint on the CONTROL instrument.
 - a. Press and hold the **set** key.
 - b. Use ▲ key and ▼ key to enter the operating temperature setpoint.
 - c. Release the set key.

WARNING:

Do not use flammable solvent or flammable material in this oven. Do not process closed containers of any substance or liquid in this oven because they may explode under heat.

Startup (Cont.)

- 3. Set HI-LIMIT instrument to a temperature 10°C to 14°C higher than the setpoint or to a temperature that should not be exceeded in the process.
- 4. Turn heater rocker switch to the ON position. The white light will also come on, indicating a heat on condition. When oven temperature reaches setpoint, the white light will cycle on and off as controlled by the CONTROL instrument.
- 5. After heating cycle is complete, turn the heater rocker switch to the OFF position. Do not turn the power off until the oven chamber temperature is below 100°C (212°F).

Optimum performance for **specific** operating conditions may require calibration. See the Calibration section in this manual.

NOTE:
If the HI-LIMIT instrument is exceeded the heater will shut down. Reset the HI-LIMIT by pushing the black button on the HI-LIMIT instrument.

Maintenance

Do not attempt any service on this oven before opening the main power disconnect switch.

Checklist

- ✓ Keep equipment clean. Gradual dirt accumulation retards air flow. A dirty oven can result in unsatisfactory operation such as unbalanced temperature in the work chamber, reduced heating capacity, reduced production, overheated components, etc. Keep the walls, floor and ceiling of the oven work chamber free of dirt and dust. Floating dust or accumulated dirt may produce unsatisfactory work results. Keep all equipment accessible. Do not permit other materials to be stored or piled against it.
- ✓ Protect controls against excessive heat. This is particularly true of controls, motors or other equipment containing electronic components. Temperatures greater than 51.5°C (125°F) should be avoided.
- ✓ Establish maintenance & checkup schedules. Do this promptly and follow the schedules faithfully. Careful operation and maintenance will be more than paid for in continuous, safe and economical operation.
- ✓ <u>Maintain equipment in good repair.</u> Make repairs immediately. Delays may be costly in added expense for labor and materials and in prolonged shut down.
- ✓ Practice safety. Make it a prime policy to know what you are doing before you do it. Make CAUTION, PATIENCE, and GOOD JUDGEMENT the safety watchwords for the operation of your oven.
- ✓ <u>Lubrication</u>, Fan motor bearings are permanently lubricated. All door latches, hinges, door operating mechanisms, bearing or wear surfaces should be lubricated to ensure easy operation.

Tests

Tests should be performed carefully and regularly. The safety of personnel as well as the condition of equipment may depend upon the proper operation of any one of the functions of these controls. Test the CONTROL instrument every 40 hours. Check that the CONTROL instrument heater LED is cycling on and off, indicating that the heater is working.

Test the HI-LIMIT instrument every 40 hours. With the oven operating at a given temperature, gradually turn the HI-LIMIT instrument knob down to the setpoint operating temperature. The HI-LIMIT instrument has tripped when the white LED above the POWER switch shuts off. Push the HI-LIMIT instrument RESET button after returning the thermostat to its original setting.

WARNINGS:

Failure to heed warnings in this manual and on the oven could result in death, personal injury or property damage.

WARNING:

Disconnect the main power switch or power cord before attempting any repairs or adjustments.

Replacement

Parts

To order or return parts, contact the Service Products Division at Despatch. The Service Products Division features our Response Center for customer service. When returning parts, the Despatch representative will provide you with an MRA (Material Return Authorization) number. The MRA number must be attached to the returned part for identification. When you are ordering parts, be sure to give the model number, serial number and the part number. This will expedite the process of obtaining a replacement part.

When you have a service need, just contact the **Response** Center at 1-800-473-7373: FAX 612-781-5353.

WARNING:

Disconnect the main power switch or power cord before attempting any repair or adjustment.

CONTROL Instrument

(Tools needed: one quarter (1/4) inch socket set screwdriver)

- 1. Disconnect power.
- 2. Remove screws with ¼ inch socket from the face of the control panel and slide it forward.
- 3. Remove wires from the old control instrument, noting which numbered wires connect to which terminals.




Figure 5 illustrates the connections to the ${\tt CONTROL}$ instrument.

- 4. Disconnect the CONTROL mounting bracket.
- 5. Remove old CONTROL instrument from control panel.
- 6. Install new CONTROL instrument into the control panel.
- 7. Secure CONTROL mounting bracket.
- 8. Reattach wires to the new CONTROL instrument. Make sure that the wires are connected correctly.
- 9. Replace control panel.

HI-LIMIT Instrument

(Tools needed: one quarter (1/4) inch socket set screwdriver)

- 1. Disconnect power.
- 2. Remove screws with ¼ inch socket from the face of the control panel and slide it forward.
- 3. Remove wires from the old HI-LIMIT instrument, noting which numbered wires connect to which terminals.




Figure 6 illustrates the connections to the HI-LIMIT instrument.

- 4. Remove mounting screws holding the HI-LIMIT instrument to the control panel.
- 5. Remove old HI-LIMIT instrument from control panel.
- 6. Install new HI-LIMIT instrument into the control panel.
- 7. Replace mounting screws.
- 8. Reattach wires to the new HI-LIMIT instrument. Make sure that the wires are connected correctly.
- 9. Replace control panel.

Heater Unit

(Tools needed: Crescent wrench, screwdriver, one quarter (1/4) inch socket set)

- 1. Remove floor plate.
 - a. Remove screws from the floor plate.
 - b. Lift floor plate out of the oven.
- 2. Disconnect heater leads from heater element with wrench. Note which wires go on which terminals.
- 3. Unscrew screws holding the heater frame to the oven body.
- 4. Remove heater and discard.
- 5. Screw down new heater frame.
- 6. Attach heater leads to appropriate terminals.
- 7. Replace interior floor and screws.

Fan Motor

(Tools needed: screwdriver, 5/32 inch Allen wrench, one quarter $(\frac{1}{4})$ inch socket set)

- 1. Remove floor plate.
 - a. Remove screws from the floor plate.
 - b. Lift floor plate out of the oven.
- 2. Loosen set screws (2) on fan wheel in middle of oven bottom. Access the fan wheel through the heater or by disconnecting and removing the heater to gain access to the fan. Refer to the Heater Unit instructions.
- 3. Place oven on its back.
- 4. Remove bottom plate. This will reveal the fan motor.
- 5. Remove fan motor.
 - a. Unscrew screws (4) holding motor mounts to body.
 - b. Disconnect motor leads from the terminal strip on power panel.
 - c. Lift the fan motor from the oven body.

After running at temperature, the fan wheel will stick to the shaft. Some force may be required to separate the fan wheel from the fan motor shaft.

- 6. Take motor mounts off old motor.
- 7. Put motor mounts onto new motor.

Fan Motor (Cont.)

- 8. Replace fan motor.
 - a. Insert shaft into shaft collar. Put fan wheel onto shaft from inside of oven.
 - b. Reattach motor mounts to oven body, making sure grommets are in place.
 - c. Reattach motor lead wires to terminal strip.

Table 6 Motor Leads to Terminal Strip

Line			Notor Leads	
Voltage	Line 1	Line 2	Tied Together*	Not Used**
120 Volt	Orange	Blue & White	Red & Black	None
240 Volt	Orange	White	Blue & Red	Black

^{*} These leads should have a wire nut to tie them together.

- 9. Replace oven bottom.
- 10. Turn oven right side up.
- 11. Adjust fan wheel for 3/16 inch clearance between wheel and inlet ring.
- 12. Tighten set screws on the fan wheel.
- 13. Check that set screws hit the flats machined into the motor shaft.
- 14. Replace floor plate.

^{**} This lead should have a wire nut on the end for isolation.

APPENDIX

Special Instructions

The oven has been tested and preset at the factory for normal operating conditions. In most applications, it will not be necessary to alter the oven's settings. This section contains additional information and reference material to access the Parameter Programming mode. This section also covers calibration for the CONTROL and for the HI-LIMIT.

The CONTROL instrument was carefully programmed at the factory using the Parameter Programming mode. The parameters that may be accessed include temperature setting, display functions and thermocouple selection.

The Calibration section covers the procedure necessary to recalibrate the CONTROL. Recalibration may be necessary if the CONTROL does not comply with known standards or to align the CONTROL for a specific operating condition. The Calibration section also covers the procedure for aligning the HI-LIMIT thermostat when HI-LIMIT recalibration is necessary.

Parameter Programming Mode

The instrument and control parameters are set through the Parameter Programming mode. In most applications, it is not necessary to alter the oven settings. The following instructions describe how to access, view and, if desired, change the parameters.

Once the Parameter Programming mode is accessed, the output LED will start blinking on and off. If a particular setting is not allowed, the display will flash. The CONTROL will not allow the display to be altered improperly.

The CONTROL will automatically exit the Parameter Programming mode if no keys are pressed for about 20 seconds. During programming, the output LED should be blinking

1. Press the **prg** key, **set** key and the **ewtr910** key simultaneously.




Figure 7 illustrates the control instrument.

- 2. Check that the LED is flashing, indicating that the control is in the program mode.
- 3. Press the A key until the desired parameter is displayed. See Program Mode Outline.

WARNING:

Make sure you understand what you are changing before doing so. Changing the program parameters will alter the functions of the CONTROL.

Parameter Programming Mode (Cont.)

- 1. Press and hold the **set** key to view the parameter setting.
- 2. While holding the **set** key, use the ▲ key and the ▼ key to change the parameter setting to the desired setting.
- 3. Press the **prg** key, **set** key and **ewtr910** key simultaneously to exit the Parameter Programming mode. The CONTROL will revert back to its normal mode.

Table 7 Parameter Program Mode Outline

Code	Name	Settings
LS1	Lower Set 1 (degrees)	50°C
HS1	Higher Set 1 (degrees)	204°C
Pb*	Proportional Band	7°C
It	Integral Time	120s
dt	Derivative Time	Os
Sr	Sampling Rate	1s
rSt	Manual Reset	0°C
Ar	Anti Reset Wind-up	20°C
od	Output Delay	0
Ct	Cycle Time	20s
drb	Dynamic Restart Band	0
dSi	Dynamic Set Increment	1
dSt	Dynamic Set Time	1
CAL	Calibration	0
Ft	Function Type	Pi
PSE	Probe Selection	FE
HC1	Heat/Cool Output	Н
rP1	Relay Protection 1	ro
LF1	LED Function 1	di
rou	Readout (degrees)	°C
dro	Display Readout	P
hdd	Half-Digit Display	n
tAb	Table of Parameters	not adjustable

^{*} If Pb is not displayed the Function Type must be first set to Pi.

Parameter Programming Mode (Cont.)

- LS1 Lower Set 1 This is the lower setpoint limit below which the user cannot change the setpoint.
- HS1 Higher Set 1 The maximum setpoint limit for chamber. The user cannot set the setpoint above the maximum setpoint.
- Pb Proportional Band Expressed in degrees. This value determines the band width on both sides of the setpoint within which the control provides proportional control.
- It Integral Time Expressed in seconds. This parameter corrects for errors in actual temperature versus the setpoint.
- Dt Derivative Time Expressed in seconds. This effect of the derivative time is in direct proportion to the time setting.
- Sr Sampling Rate Expressed in seconds. Time between two successive read-outs. Recommended setting is 1.
- rSt Manual Reset Expressed in degrees. This allows the proportional band to be moved up or down.
- Ar Anti-reset Windup This is the half-band in which the integral action takes place. The higher the setting, the stronger the integral action.
- od Output Delay This provides a delay selection for the outputs in applications where noise may cause brief erroneous signals from the sensor to the controller.
- Ct Cycle Time Expressed in seconds. This is the total time on one ON/OFF cycle of the relay during the proportional action.
- drb Dynamic Restart Band This is a soft start function. When temperature falls below or rises above this restart band, another soft cycle is initiated. The value of this parameter represents half of the total band.
- dSi Dynamic Set Increment Expressed in degrees of the setpoint. This parameter represents the dynamic increase of the setpoint. See Dynamic Restart Band. A setting of 0 disables this function.

Parameter Programming Mode (Cont.)

- dSt Dynamic Set Time Time value between two successive dynamic setpoint increases.
- CAL Calibration The number of degrees the control will offset the display from the sensor input. This parameter is used to align the oven's actual chamber temperature with the display appearing on the control. The ± sign determines whether the adjustment is made upward or downward. See Calibration at the end of this addendum.
- Ft Function Type Control mode selection.
 on = ON/OFF control,
 Pi = PID control
 This oven operates only with the Pi control setting.
- PSE Probe Selection Input type. FE = J T/C, Cr = K T/C, rh = S T/C
- HC1 Heat/Cool Output Heating = H, Cooling = C. Set to Heating (H) for heating applications.
- rP1 Relay Protection 1 Determines the status of the relay in the event of a sensor error.

 ro = relay open (factory setting)

 rc = relay closed

 Use the ro setting for heating applications.
- LF1 LED Function 1 Determines whether the light is on or off when the relay is energized.

 di = direct, light is ON when output relay is energized in = reverse, light is OFF when the output relay is energized
- rou Readout (degrees)
 C = degrees celsius
 F = degrees fahrenheit
- dro Display Readout P = Process value S = Setpoint value
- hdd Half-Digit Display The right most digit can be set to read-out in 0 or 5 only automatically rounding to the nearest value.
- tAb Table of Parameters Factory setting, cannot be changed.

Calibration

CONTROL

The CONTROL instrument has been tested and calibrated at the factory. Under normal operating conditions, recalibration should not be necessary. However, if the instrument does not comply with known standards, OR if the user would like to recalibrate the CONTROL for a specific operating condition, then recalibration is easily accomplished.

Calibration Instructions

(Equipment needed: Temperature Measuring Device with a Compatible Temperature Sensor)

- Verify that the CAL (Calibration) programmed in the CONTROL is 0. Refer to Instructions on viewing the parameter in the OPERATING mode in the Appendix of this manual
- 2. Locate the temperature sensor of the temperature measuring device at the center of the chamber.
- 3. Operate the oven until it reaches the desired operating temperature and the CONTROL is regulating. The user may wish to have a loaded chamber with a standard amount of product to simulate a specific operating condition. It will take several minutes for the unit to stabilize at the controlled temperature. Allow at least 30 minutes of operation at the stabilized temperature before proceeding.
- 4. Subtract the average controlled temperature (number appearing on the CONTROL display) from the actual oven temperature (number appearing on the temperature measuring device display). The CONTROL and the device must be in the same scale (°C or °F).

Actual Oven Temperature - Controlled Temperature = calculated value

5. Enter the calculated value from Step 4 as the new CAL (Calibration) value in the instruments.

HI-LIMIT

The HI-LIMIT equipped on the ovens is a non-indicating HI-LIMIT (a device that does not display sensor temperature). The HI-LIMIT has been tested and calibrated at the factory and recalibration should not be necessary. However, since the HI-LIMIT is a non-indicating device, the HI-LIMIT may need to be recalibrated against known standards so that the temperature dial indication on the HI-LIMIT matches the HI-LIMIT thermocouple input.

Calibration Instructions

(Tools Needed: 1/16 inch Allen wrench)

- 1. Set CONTROL instrument at 14°C (25°F) above the desired operating temperature.
- 2. Push black button to reset the HI-LIMIT instrument and operate oven until the CONTROL instrument is regulating. Allow the oven to soak at temperature for 30 minutes.
- 3. Adjust the HI-LIMIT instrument downward until it trips (HEATER ON LED goes off).
- 4. Loosen the set screw on the black knob on the HI-LIMIT with a 1/16 inch Allen wrench.
- 5. Align the HI-LIMIT dial with the CONTROL temperature display.
- Tighten the set screw.
- Reset CONTROL instrument at the desired operating temperature. The two instruments are now set in their correct positions.

It will be necessary to reset the HI-LIMIT instrument whenever it has tripped. The HI-LIMIT instrument may be reset by first allowing the oven chamber to cool slightly (or turning the HI-LIMIT instrument thermostat up several degrees) and pushing the black reset button.

Accessories

The LDB Series Ovens have options that can easily be field installed.

Table 8 Accessories

Option	Functional Description		
Recorder kit	The round chart recorder follows the temperature changes and records them for permanent record.		
240 Volt conversion kit	Converts 120 volt oven for 240 volt operation. Available on LDB 1-17, LDB 1-24, LDB 1-43 and LDB 1-69. 240 V conversion must be factory installed to conform to UL listing requirements.		
Process timer kit	Shuts down the heater at the end of the heating cycle. Available in 6 or 12 hour ranges.		
Running time meter kit	Logs the overall processing time. The digital meter charts up to 99,999.9 hours of process time and cannot be reset.		
Signal timer kit	One hour timer signals the end of the cycle with a bell.		
Extra shelves			
Base cabinets			
Stacking kit			

The above items can all be field installed. For further information on these items or other available options, please contact your Despatch representative.

Troubleshooting

Equipment which operates for long periods of time may develop problems. Below are possible problems and suggested solutions. If you have a problem not listed and do not know what to do, contact Despatch Industries at our toll free Help Line 800-473-7373.

Difficulty	Probable Cause	Suggested Remedy
Failure to heat	No power	Check power source and/or oven and wall fuses.
	Broken or frayed cord	Replace with new cord and plugset.
	Burned out heater(s)	Replace element (see warranty statement).
	CONTROL instrument malfunction	Replace controller.
	Loose wire connections	Disconnect power and check connections behind control panel.
Slow heat up	Improperly loaded workload	Reduce load or redistribute load in chamber.
	Low line voltage	Supply sufficient power and proper connections. Check for circuit overload.
	1 or 2 heating elements burned out	Check heater amperage on the nameplate. Replace burned out element (see warranty statement, back page).
	Vent(s) are wide open.	Close vent(s).
	Fan motor failure.	Replace fan motor.
Frequent heater element out	Harmful fumes generated by load	Increase vent opening or discontinue process.
	Spillage or splattering of material on heater elements	Disconnect power and clean oven chamber and elements.
	Overheating oven	Check the HI-LIMIT instrument.
Erratic temp.	CONTROL instrument malfunction	Check control parameters before replacing the CONTROL instrument.

Troubleshooting (Cont.)

Difficulty	Probable Cause	Suggested Remedy
Inaccurate temp.	CONTROL instrument miscalibration	Recalibrate CONTROL instrument.
Excess surface or door temp.	Door seal deterioration	Replace door seal.
Improper airflow	Fan motor failure	Replace fan motor.
	Unbalanced fan wheel	Replace fan wheel.
Excessive	Dirty fan wheel	Clean fan.
vibration	Unbalanced fan wheel	Replace fan wheel.
Grinding noise	Fan wheel has shifted or fallen.	Inspect the wheel. Wheel should be 3/16 of an inch lower than the top of the fan housing.
Oven will not control at	HI-LIMIT instrument set too low	Set the HI-LIMIT higher.
setpoint	CONTROL instrument malfunction	Check control parameters before replacing CONTROL instrument.
	Air friction of recirculation fan	Open exhaust air vent. Unit will not control at minimum operating temperature with vent(s) closed.
Heater does not shutdown until temp. reaches the	CONTROL instrument malfunction	Check voltage across the input terminals of the SSR relay. Replace CONTROL instrument if 5VDC output exists.
HI-LIMIT setting	SSR Relay malfunction	Replace relay if no output exists.

Warranty

For years Despatch has delivered an exceptional product backed by a strong sense of responsibility and drive for long term customer satisfaction. These business principles enable us to offer the exclusive and comprehensive "Classic™ Four Plus One Service Warranty Program".

Despatch Classic™ Service Warranty Program

The basis of this outstanding, exclusive service program is a four-year replacement warranty covering defects in workmanship or material on all Despatch manufactured components and assemblies. Add to this a comprehensive 1 year parts warranty on the entire unit and you have one of the strongest warranties in the industry.

Immediate Service Response

The key to the Classic™ Service Program is response. A toll free Help-Line connects you to our Customer Service response center giving you immediate access to specialized assistance. Our customer service Product Service Technicians have over 200 years experience and access to detail design and manufacturing documentation specific to your Despatch unit. This exacting level of service is a benefit only Despatch can provide and means that you can expect speedy, accurate and the most cost effective response.

Field Service Network

A growing network of Service Professionals are available to support your Despatch equipment. From routine calibration and preventive maintenance to emergency breakdown response, our service network is positioned to reach 90 percent of our installed base within four hours. This is service you can depend on.

Best Service Protection in the Industry!




Figure 8 illustrates the Classic [™]4 Plus 1 Service Program

^{**} Receive 5th year parts replacement warranty free with purchase of PM+ or Total Service Agreement within the first two years.

Classic™ Series Warranty

Parts and Material

Despatch warrants all parts and assemblies manufactured by Despatch for the ClassicTM Series oven to be free from defects in material and workmanship for a period of four (4) years from the date of shipment or start-up, by Authorized Despatch Service Representative, whichever is later.

Despatch further warrants all parts and assemblies to be free from defects in material and workmanship for a period of one (1) year from date of shipment or start-up, by Authorized Despatch Service Representative, whichever is later.

Despatch will repair or replace, at our option, f.o.b. Despatch's factory, parts covered by this warranty. Despatch is not responsible for parts defects resulting from misuse, abuse, acts of nature or utility performance not to Despatch specification including electrical, environmental and fresh air/exhaust provisions.

Labor and Expenses

Despatch ClassicTM Series warranty cover parts replacement or repair. Labor and other expenses related to the removal and replacement of such parts are the owners responsibility as is any necessary reprogramming, calibration and certification.

Exclusions/Limitations of Liability

The foregoing warranty shall be deemed valid and binding upon Seller if and only if Purchaser loads, operates and maintains the equipment supplied hereunder in accordance with the instruction manual provided upon delivery of the equipment. Seller does not guarantee the process of manufacture by Purchaser or the quality of product to be produced by subject equipment. This warranty does not cover expenses to diagnose, repair or replace components or associated failures.

Exclusions/Limitations of Liability (Cont.)

Parts failures caused by improper operation, abuse, misuse, acts of nature, and nonconforming utilities and environments are not covered by this warranty.

Despatch shall not in any event be liable for indirect, special, consequential or liquidated damages or penalties, including loss of revenue, profits or business opportunities resulting from interruption of product production. Despatch shall further be held blameless for any damages or expenses resulting from delays in our attempts to diagnose and repair the equipment, unavailability of spare parts or inaccessibility of the equipment. Specifically excluded from this warranty is responsibility for internal and external corrosion damage to the equipment.

Emergency Service

In an emergency situation, Customer agrees to:

- 1. Immediately shut off fuel or energy supply (gas and electricity).
- 2. Call 911 for emergency assistance if needed.
- 3. Call Despatch Service Help-Line at 800-473-7373.

Non-Compliance

Despatch reserves the right to suspend and withhold service as provided under this Warranty in the event of non-compliance by the Customer to any terms and conditions of this Warranty. Further, Despatch is held harmless for any loss of production, incurred expenses, or other inconveniences due to suspension of service under this non-compliance provision.

THE FOREGOING WARRANTY IS EXCLUSIVE AND IN LIEU OF ALL OTHER EXPRESS AND IMPLIED WARRANTIES WHATSOEVER, AND SPECIFICALLY THERE ARE NO IMPLIED WARRANTIES OF MERCHANTABILITY OR OF FITNESS FOR A PARTICULAR PURPOSE.

SEE ATTACHED DESCRIPTION OF DESPATCH CUSTOMER SERVICE PROCEDURES AND RELATED CUSTOMER RESPONSIBILITIES WHICH ARE INTEGRAL TO THIS WARRANTY.

THE FOREGOING WARRANTY IS NOT TRANSFERRABLE IN SITUATIONS WHERE EQUIPMENT OWNERSHIP IS TRANSFERRED TO ANOTHER PARTY.

Despatch Customer Service

Procedures and Customer Responsibilities

To provide the most effective service to our customers under this warranty, all requests for repairs are to be initiated by the Customer by telephone to the Despatch Service Help Line, 800-473-7373. The Standard Period of Maintenance (SPM) is defined as 8 a.m. to 5 p.m. local time, excluding weekends and Despatch Holidays. Calls placed within the SPM will be handled as follows.

Help Line calls connect the customer with the Despatch Response Center. The Response Center will record all pertinent information, including SERIAL and MODEL NUMBER of the unit(s), the urgency and nature of the problem, and the name and phone number of the caller or other contact. This information will be passed to the first available service support technician who will research the units serial file so as to be familiar with customer unit when he calls the customer back. Despatch service technicians will make every effort to call back within four (4) working hours, or less, from receipt of the initial call. Despatch will advise the Customer on suggested steps and/or tests to either resolve the problem or help to confirm the diagnosis of the problem.

Customer Agrees to cooperate in performing such tests and attempting to resolve the problem as quickly as possible. Customer also agrees to replace minor parts such as fuses, latches, etc. as instructed by Despatch Service Technicians. This approach has Despatch and the Customer cooperating to effect the most expedient and cost effective repair and minimize down time. If in Despatch's sole judgement, the equipment cannot be repaired in this manner, an on-site visit by a Despatch authorized service representative may be scheduled to repair the equipment. Customer agrees that, when requested and authorized, such charges will be paid by the Customer within 30 days from receipt of invoice.


Attachment A - Sustained Service Support

At Despatch long term customer satisfaction means more than responding quickly and effectively to our customers' service requirements. It means offering comprehensive customer support well beyond the scope of our initial contractual commitment. Despatch's Service Products Division offers a Total Service Agreement package or a Preventive Maintenance Plus agreement (PM+). These service products are unique in the industry and offer the following benefits to our customer:


- Priority response for minimum production interruption.
- Preventive maintenance for longer product life.
- Discount on parts and services where applicable.
- Single payment for reduced billing expense.
- Elimination of need for a separate purchase order for each service requested.

Because these extra service options are aimed at extending our new equipment productivity, we will also extend the Despatch four year manufactured parts warranty for another 12 months. This bonus warranty is automatically yours when you purchase a service agreement from Despatch within the first 12 months after shipment of the equipment.

Drawings



LDB1-17 240V




Item	Part Number	Quantity	Description
CL	098849	1	Despatch CONTROL
HL	086420	1	Despatch HI-LIMIT
1F	007453	1*	250 volt 6 amp fuse
M	008333	1	1/25 HP motor
С	031233 (015115)	1	Power cord 120V (240V)
1,2SW	074116	2	Switch gray DPST
1LT	008657	1	Amber pilot light
2LT	019150	1	White pilot light
1CR	080187	1	20 amp relay
2CR	073739	1	100-A09ND3 contactor
Н	007818	1	1200 watt heater
1,2T/C	111322	2	T/C LEB/LDB-4 72 inches long test
Т	012479	1	50VA transformer
QARC	046784	1	104M06QC100 quencharc
VAR	014643	2	V130LA10A Varistor
CF	006049 (015229)	1	MU2A1 fan 120VAC (MU3A1 fan 240VAC)


^{*} Two fuses required on 240V

LDB1-24 120V

120V-1PH-60HZ




LDB1-24 240V




Item	Part Number	Quantity	Description
CL	098849	1	Despatch CONTROL
HL	086420	1	Despatch HI-LIMIT
1F	007453	1*	250 volt 6 amp fuse
M	008333	1	1/25 HP motor
С	031233 (015115)	1	Power cord 120V (240V)
1,2SW	074116	2	Switch gray DPST
1LT	008657	1	Amber pilot light
2LT	019150	1	White pilot light
1CR	080187	1	20 amp relay
2CR	073739	1	100-A09ND3 contactor
Н	007818	1	1200 watt heater
1,2T/C	111322	2	T/C LEB/LDB-4 72 inches long test
T	012479	1	50VA transformer
QARC	046784	1	104M06QC100 quencharc
VAR	014643	2	V130LA10A Varistor
CF	006049 (015229)	1	MU2A1 fan 120VAC (MU3A1 fan 240VAC)

^{*} Two fuses required on 240V

LDB1-43 120V



LDB1-43 240V




Connection Legend


Item	Part Number	Quantity	Description
CL	098849	1	Despatch CONTROL
HL	086420	1	Despatch HI-LIMIT
1F	007453	1*	250 volt 6 amp fuse
M	008333	1	1/25 HP motor
С	074116	1	Power cord
1,2SW	074116	2	Switch gray DPST
1LT	008657	1	Amber pilot light
2LT	019150	1	White pilot light
1CR	080187	1	20 amp relay
2CR	073739	1	100-A09ND3 contactor
Н	054353	1	1600 watt heater
1,2T/C	111322	2	T/C LEB/LDB-4 72 inches long test
QARC	046784	1	104M06QC100 quencharc
VAR	014643	2	V130LA10A Varistor
CF	006049	1	MU2A1 fan 120VAC

^{*} Two fuses required on 240V

LDB1-69 120V



LDB1-69 240V




Item	Part Number	Quantity	Description
CL	098849	1	Despatch CONTROL
HL	086420	1	Despatch HI-LIMIT
1F	007617	1*	250 volt 30 amp fuse
2F	007453	1	250 volt 6 amp fuse
М	008333	1	1/25 HP motor
1,2SW	074116	2	Switch gray DPST
1LT	008657	1	Amber pilot light
2LT	019150	1	White pilot light
1CR	080187	1	20 amp relay
2CR	073739	1	100-A09ND3 contactor
Н	007818	3	1200 watt heater
1,2T/C	111322	2	T/C LEB/LDB-4 72 inches long test
Т	012479	1	50VA transformer
QARC	046784	1	104M06QC100 quencharc
VAR	014643	2	V130LA10A Varistor
CF	006049 (015229)	1	MU2A1 fan 120VAC (MU3A1 fan 240VAC)**

^{*} Two fuses required on 240V


** Depends with voltage

LDB2-18 240V



Item	Part Number	Quantity	Description
CL	098849	1	Despatch CONTROL
HL	086420	1	Despatch HI-LIMIT
1F	007615	2	250 volt 25 amp fuse
2F	007453	2	250 volt 6 amp fuse
M	008333	2	1/25 HP motor
1,2SW	074116	2	Switch gray DPST
1LT	008657	1	Amber pilot light
2LT	019150	1	White pilot light
1CR	080187	1	20 amp relay
2CR	073739	1	100-A09ND3 contactor
Н	007818	2	1200 watt heater
1,2T/C	111322	2	T/C LEB/LDB-4 72 inches long test
Т	012479	1	50VA transformer
QARC	046784	1	104M06QC100 quencharc
VAR	014643	2	V130LA10A Varistor
CF	015229	1	MU3A1 fan 120VAC

LDB2-27 240V



Item	Part Number	Quantity	Description
CL	098849	1	Despatch CONTROL
HL	086420	1	Despatch HI-LIMIT
1F	007617	2	250 volt 30 amp fuse
2F	007453	2	250 volt 6 amp fuse
М	008333	2	1/25 HP motor
1,2SW	074116	2	Switch gray DPST
1LT	008657	1	Amber pilot light
2LT	019150	1	White pilot light
1CR	080187	1	20 amp relay
2CR	073739	1	100-A09ND3 contactor
Н	007819	2	2400 watt heater
1,2T/C	111322	2	T/C LEB/LDB-4 72 inches long test
Т	012479	1	50VA transformer
QARC	046784	1	104M06QC100 quencharc
VAR	014643	2	V130LA10A Varistor
CF	015229	1	MU3A1 fan 120VAC